

Exercise #21. SIM Agar. Determine H₂S, indole production and motility. p. 54

| | | |
|---|---|--|
| <p>Day 1 Mon 7/20/09</p> | <p>Procedure:</p> <ul style="list-style-type: none"> • work in groups of 4 • stab inoculate 7 SIM agar tubes <ul style="list-style-type: none"> ○ 3 known controls ○ 4 unknowns (1 for each person in your group) • INCUBATE: 37°C; 48 hours | <p>Organisms: Controls:</p> <ol style="list-style-type: none"> 1. <i>E. coli</i> 2. <i>P. vulgaris</i> 3. <i>K. pneumoniae</i> <p>Your unknowns:</p> |
| <p>Day 2 Wed 7/22/09</p> | <p>Results:</p> <ul style="list-style-type: none"> • observe tubes • H₂S <ul style="list-style-type: none"> ○ presence of black precipitate • Motility, examine stab line <ul style="list-style-type: none"> ○ if not motile, cells will not leave line ○ if motile, cells will move off line • Indole <ul style="list-style-type: none"> ○ add Kovacs reagent to each tube ○ if indole present, it will move up into the Kovacs reagent ○ cherry red color will form • record your results on p. 31-2. 55 • answer questions • refer to dichotomous key on p 33-34 | <p>NOTE:</p> <p>Record H₂S and motility BEFORE testing for indole</p> |



